


Anti-Nlrp12 Polyclonal Antibody For Wb, Elisa - P59046

Item Number: CM0014326



Introduction

Anti-NLRP12 rabbit polyclonal antibody validated for WB and ELISA in human, with rat cross-reactivity. NLRP12 is a key inflammation mitigator, inhibiting NF-κB, NOD2, and RIGI signaling. Ideal for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	NLRP12 Rabbit pAb
Remarks/Alias	RNO; PAN6; RNO2; FCAS2; NALP12; PYPAF7; CLR19.3; NLRP12
Species	Human
Gene ID	91662
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 862-1061 of human NLRP12 (NP_653288.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P59046
Protein Weight	120kDa
Shipping	Ice bag