


Anti-Nlrp12 Rabbit Polyclonal Antibody For Wb And Elisa - P59046

Item Number: CM0029140



Introduction

High-quality rabbit polyclonal antibody against human NLRP12 (Monarch-1, P59046) validated for WB and ELISA. Ideal for studying inflammation, NF-κB signaling, and immune regulation in leukocytes.

[Learn More](#)

Parameter	Specification
Product Name	NLRP12 Rabbit pAb
Remarks/Alias	RNO; PAN6; RNO2; FCAS2; NALP12; PYPAF7; CLR19.3; NLRP12
Species	Human
GeneID (Human)	91662
GeneID	91662
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P59046
Protein Weight	120kDa
Shipping	Ice bag