

Anti-Nlrp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nx02

Item Number: CM0028639



Introduction

Rabbit polyclonal antibody against human NLRP2, suitable for WB, IF/ICC, and ELISA. Detects endogenous NLRP2 in human and mouse samples. Useful for inflammasome and NF-kB signaling research.

[Learn More](#)

Parameter	Value
Product Name	NLRP2 Rabbit pAb
Remarks/Alias	NBS1; PAN1; NALP2; PYPAF2; CLR19.9; NLRP2
Species	Human
GeneID (human)	55655
GeneID	55655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 665-820 of human NLRP2 (NP_060322.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NX02
Protein Weight	121kDa
Shipping	Ice bag