


Anti-Nlrp3 Rabbit Monoclonal Antibody For Wb, Elisa - Q96P20

Item Number: CM0010196



Introduction

Rabbit monoclonal antibody targeting human NLRP3, a key inflammasome sensor. Validated for Western blot and ELISA. Reacts with human samples. Molecular weight ~118 kDa.

[Learn More](#)

Parameter	Value
Product Name	NLRP3 Rabbit mAb
Remarks/Alias	Ali; AVP; FCU; MWS; FCAS; KEFH; CIAS1; FCAS1; NALP3; C1orf7; CLR1.1; DFNA34; PYPAF1; AGTAVPRL
Species	Human
GeneID (Human)	114548
GeneID	114548
Immunogen	Recombinant protein of human NLRP3
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96P20
Protein Weight	118kDa
Shipping	Ice bag