


Anti-Nlrp3 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q96P20

Item Number: CM0012758



Introduction

Rabbit monoclonal antibody against human NLRP3, validated for Western blot, immunoprecipitation, and ELISA. Essential tool for studying NLRP3 inflammasome activation, IL-1 β /IL-18 secretion, pyroptosis, and Th2 differentiation.

[Learn More](#)

Parameter	Value
Product Name	NLRP3 Rabbit mAb
Remarks/Alias	Ali; AVP; FCU; MWS; FCAS; KEFH; CIAS1; FCAS1; NALP3; C1orf7; CLR1.1; DFNA34; PYPAF1; AGTAVPRL
Species	Human
Gene ID (Human)	114548
Gene ID	114548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-210 of human NLRP3 (NP_004886.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q96P20
Protein Weight	118 kDa
Shipping	Ice bag