


Anti-Nlrp3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96P20

Item Number: CM0015043



Introduction

High-quality rabbit polyclonal antibody targeting NLRP3 (cryopyrin), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for inflammasome activation and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	NLRP3 Rabbit pAb
Remarks/Alias	AIL; AVP; FCU; MWS; FCAS; KEFH; CIAS1; FCAS1; NALP3; C1orf7; CLR1.1; DFNA34; PYPAF1; AGTAVPRL; NLRP3
Species	Human
GeneID (human)	114548
GeneID	114548
Immunogen	Recombinant Protein A synthetic peptide corresponding to a sequence within amino acids 1-220 of human NLRP3 (NP_001120934.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P20
Protein Weight	118kDa
Shipping	Ice bag