

Anti-Nlrp5 Polyclonal Antibody For Wb, Elisa - P59047

Item Number: CM0025773



Introduction

Rabbit polyclonal antibody targeting NLRP5 (MATER, NALP5), validated for Western blot and ELISA applications in human, mouse, and rat samples. Suitable for studying subcortical maternal complex and early embryonic development.

[Learn More](#)

Parameter	Value
Product Name	NLRP5 Rabbit pAb
Remarks/Alias	MATER, NALP5, PAN11, PYPAF8, CLR19.8, NLRP5
Species	Human
GeneID	126206
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 157-256 of human NLRP5 (NP_703148.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P59047
Protein Weight	134kDa
Shipping	Ice bag