


Anti-Nlrp6 Rabbit Monoclonal Antibody For Wb, Elisa - P59044

Item Number: CM0010398



Introduction

High-quality rabbit monoclonal antibody targeting human NLRP6, validated for WB and ELISA applications. Suitable for research on inflammasome activation, innate immunity, and host defense mechanisms. Swiss-Prot: P59044.

[Learn More](#)

Parameter	Value
Product Name	NLRP6 Rabbit mAb
Remarks/Alias	AVR; NAVR; PAN3; NALP6; PYPAF5; CLR11.4; NAVR/AVR
Species	Human
Gene ID (Human)	171389
Gene ID	171389
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 793-892 of human NLRP6 (NP_612202.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P59044
Protein Weight	99 kDa
Shipping	Ice bag