

Anti-Nlrp6 Rabbit Polyclonal Antibody For Wb, Elisa - Q91Ws2

Item Number: CM0019991

Introduction

Rabbit polyclonal antibody targeting mouse Nlrp6 for WB and ELISA. Recognizes mouse and rat Nlrp6. Useful for inflammasome biology, gut mucosal immunity, and host-pathogen interaction studies.

[Learn More](#)

Parameter	Value
Product Name	Nlrp6 Rabbit pAb
Remarks/Alias	Avr; Navr; Nalp6; Pypaf5; Non-AVR; Navr/Avr; Nlrp6
Species	Mouse
GeneID	101613
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 124-237 of mouse Nlrp6(NP_598707.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q91WS2
Protein Weight	97kDa
Shipping	Ice bag