


Anti-Nlrp1 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - Q86Ut6

Item Number: CM0022669



Introduction

KO validated rabbit polyclonal antibody targeting NLRX1, suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying antiviral innate immunity, autophagy, and NLRP3 inflammasome pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NLRX1 Rabbit pAb
Remarks / Alias	NOD5; NOD9; NOD26; DLNB26; CLR11.3; X1
Species	Human
Gene ID	79671
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 87-360 of human NLRX1 (NP_001269073.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UT6
Protein Weight	108kDa
Shipping	Ice bag