

Anti-Nmdar1 Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - Q05586

Item Number: CM0030132



Introduction

Rabbit polyclonal antibody targeting human NMDAR1 (GRIN1, NR1). Suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects the 105 kDa NMDA receptor subunit critical for synaptic plasticity and learning.

[Learn More](#)

Parameter	Value
Product Name	NMDAR1 Rabbit pAb
Remarks/Alias	NR1; MRD8; GluN1; NMDA1; DEE101; NDHMSD; NDHMSR; NMD-R1; NMDAR1
Species	Human
GeneID (Human)	2902
GeneID	2902
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05586
Protein Weight	105kDa
Shipping	Ice bag