


Anti-Nmdar2A Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q12879

Item Number: CM0013316



Introduction

Anti-NMDAR2A rabbit polyclonal antibody for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Targets glutamate receptor NMDA 2A (GRIN2A, 165kDa, UniProt Q12879). Ideal for synaptic plasticity studies.

[Learn More](#)

Parameter	Value
Product Name	NMDAR2A Rabbit pAb
Remarks/Alias	LKS; EPND; FESD; NR2A; GluN2A; NMDAR2A
Species	Human
Gene ID (Human)	2903
Gene ID	2903
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1130-1400 of human NMDAR2A (NP_000824.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12879
Protein Weight	165kDa
Shipping	Ice bag