


Anti-Nmdar2A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q12879

Item Number: CM0018552



Introduction

Rabbit polyclonal antibody targeting human NMDAR2A (NR2A, GluN2A) for WB, IF/ICC, ELISA applications. Cross-reacts with mouse and rat. Suitable for NMDA receptor studies in neuroscience.

[Learn More](#)

Parameter	Value
Product Name	NMDAR2A Rabbit pAb
Remarks/Alias	LKS; EPND; FESD; NR2A; GluN2A; NMDAR2A
Species	Human
Gene ID	2903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1130-1400 of human NMDAR2A (NP_000824.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q12879
Protein Weight	165kDa
Shipping	Ice bag