

Anti-Nmdar2B Rabbit Monoclonal Antibody For Wb, Elisa - Q13224

Item Number: CM0005599



Introduction

Recombinant anti-human NMDAR2B (GRIN2B) rabbit monoclonal antibody validated for Western blot (WB) and ELISA, cross-reacting with mouse and rat. Ideal for neuroscience research on NMDA receptors and synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	NMDAR2B Rabbit mAb
Remarks/Alias	NR3; MRD6; NR2B; hNR3; DEE27; EIEE27; GluN2B; NMDAR2B
Species	Human
GeneID (Human)	2904
GeneID	2904
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1133-1484 of human GRIN2B (NP_000825.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13224
Protein Weight	166kDa
Shipping	Ice bag