


Anti-Nme1/Nm23A Polyclonal Antibody For Wb, If/Icc, Elisa - P15531

Item Number: CM0020965



Introduction

Rabbit polyclonal antibody targeting human NME1/NM23A (nucleoside diphosphate kinase A). Validated for WB, IF/ICC, ELISA. Detects a 17 kDa band; cross-reacts with human and rat. Ideal for metastasis, nucleotide metabolism, and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	NME1/NM23A Rabbit pAb
Remarks/Alias	NB; AWD; NBS; GAAD; NDKA; NM23; NDPKA; NDPK-A; NM23-H1; NME1/NM23A
Species	Human
Gene ID (Human)	4830
Gene ID	4830
Immunogen	Synthetic peptide corresponding to amino acids 10-110 of human NME1/NM23A (NP_000260.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P15531
Protein Weight	17kDa
Shipping	Ice bag