


Anti-Nme1/Nm23A Rabbit Monoclonal Antibody For Wb And Elisa - P15531

Item Number: CM0011258



Introduction

Rabbit monoclonal antibody against human NME1/NM23A (NDKA, NM23-H1). Validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for research on tumor metastasis, nucleoside metabolism, and cell motility.

[Learn More](#)

Parameter	Value
Product Name	NME1/NM23A Rabbit mAb
Remarks/Alias	NB; AWD; NBS; GAAD; NDKA; NM23; NDPKA; NDPK-A; NM23-H1; NME1/NM23A
Species	Human
Gene ID	4830
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NME1/NME1/NM23A (P15531).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15531
Protein Weight	17 kDa
Shipping	Ice bag