


Anti-Nme1/Nm23A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P15531

Item Number: CM0023268



Introduction

Rabbit polyclonal antibody against human NME1/NM23A, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Targets nucleoside diphosphate kinase A, a key regulator of nucleotide homeostasis, metastasis, and DNA fragmentation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NME1/NM23A Rabbit pAb
Remarks/Alias	NB; AWD; NBS; GAAD; NDKA; NM23; NDPKA; NDPK-A; NM23-H1; 3A
Species	Human
Gene ID (Human)	4830
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 10-110 of human NME1/NM23A (NP_000260.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15531
Protein Weight	17kDa
Shipping	Ice bag