


Anti-Nme1/Nm23A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P15531

Item Number: CM0015734



Introduction

Rabbit polyclonal antibody against human NME1/NM23A (NDPK A), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Target involved in nucleotide metabolism and metastasis inhibition.

[Learn More](#)

Parameter	Value
Product Name	NME1/NM23A Rabbit pAb
Remarks/Alias	NB; AWD; NBS; GAAD; NDKA; NM23; NDPKA; NDPK-A; NM23-H1; NME1/NM23A
Species	Human
Gene ID (Human)	4830
Gene ID	4830
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15531
Protein Weight	17kDa
Shipping	Ice bag