


Anti-Nme2/Nm23B Monoclonal Antibody For Wb, Elisa - P22392

Item Number: CM0010451



Introduction

Recombinant rabbit monoclonal antibody against human NME2/NM23B (P22392). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying nucleotide metabolism, transcriptional regulation, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	NME2/NM23B Rabbit mAb
Remarks/Alias	PUF; NDKB; NDPKB; NM23B; NDPK-B; NM23-H2; NME2/NM23B
Species	Human
GeneID (Human)	4831
GeneID	4831
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P22392
Protein Weight	17kDa
Shipping	Ice bag