

Anti-Nme6/Nm23-H6 Polyclonal Antibody For Wb, Ihc-P, Elisa - O75414

Item Number: CM0013683

Introduction

Rabbit polyclonal antibody targeting NME6 (NM23-H6, IPIA-alpha). Validated for WB, IHC-P, ELISA in human, mouse, rat. Detects 21 kDa mitochondrial NDP kinase essential for mtDNA maintenance. Ideal for mitochondrial research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	NME6 Rabbit pAb
Remarks/Alias	NDK 6; NM23-H6; IPIA-ALPHA; NME6
Species	Human
GeneID (human)	10201
GeneID	10201
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 9-194 of human NME6 (NP_005784.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75414
Protein Weight	21kDa
Shipping	Ice bag