

Anti-Nme8 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N427

Item Number: CM0020076

Introduction

Rabbit polyclonal antibody targeting human NME8 (NM23-H8/TXNDC3/Sptrx-2), validated for WB and ELISA. Detects human and mouse NME8 (~67 kDa), ideal for spermatogenesis and ciliopathy research.

[Learn More](#)

Parameter	Value
Product Name	NME8 Rabbit pAb
Remarks/Alias	CILD6; DNAI8; SPTRX2; TXNDC3; NM23-H8; sptrx-2; HEL-S-99; NME8
Species	Human
GeneID (Human)	51314
GeneID	51314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 299-588 of human NME8 (NP_057700.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N427
Protein Weight	67kDa
Shipping	Ice bag