

Anti-Nmi Rabbit Monoclonal Antibody For Wb, Ihc-P And Elisa - Q13287

Item Number: CM0007427

Introduction

Rabbit monoclonal antibody against human NMI (N-myc-interactor, 35kDa), suitable for WB, IHC-P, and ELISA. Ideal for investigating innate immune signaling, STAT pathways, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	NMI Rabbit mAb
Remarks/Alias	NMI
Species	Human
GenelD (Human)	9111
GenelD	9111
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13287
Protein Weight	35kDa
Shipping	Ice bag