


Anti-Nmnat1 Polyclonal Antibody For Wb And Elisa - Q9Han9

Item Number: CM0028842



Introduction

High-affinity rabbit polyclonal antibody targeting human NMNAT1. Validated for Western blot and ELISA applications. A key tool for studying NAD⁺ metabolism and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	NMNAT1 Rabbit pAb
Remarks/Alias	LCA9, NMNAT, PNAT1, SHILCA, NMNAT1
Species	Human
GeneID (Human)	64802
GeneID	64802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human NMNAT1 (NP_073624.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HAN9
Protein Weight	32kDa
Shipping	Ice bag