


Anti-Nmnat1 Polyclonal Antibody For Wb, Elisa - Q9Han9

Item Number: CM0026921



Introduction

High-quality rabbit polyclonal antibody against human NMNAT1 (SWISS: Q9HAN9). Validated for WB and ELISA, suitable for research on NAD⁺ biosynthesis, chromatin remodeling, and neural protection.

[Learn More](#)

Parameter	Value
Product Name	NMNAT1 Rabbit pAb
Remarks/Alias	LCA9; NMNAT; PNAT1; SHILCA; NMNAT1
Species	Human
GeneID (Human)	64802
GeneID	64802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human NMNAT1 (NP_073624.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HAN9
Protein Weight	32kDa
Shipping	Ice bag