

Anti-Nmnat2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bzq4

Item Number: CM0015600

Introduction

High-quality rabbit polyclonal antibody against human NMNAT2, validated for WB and ELISA. Detects endogenous NMNAT2 in human, mouse, and rat samples. Ideal for studying NAD+ metabolism and axon maintenance.

[Learn More](#)

Parameter	Value
Product Name	NMNAT2 Rabbit pAb
Remarks/Alias	PNAT2; C1orf15; NMNAT2
Species	Human
GeneID (Human)	23057
GeneID	23057
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human NMNAT2 (NP_733820.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZQ4
Protein Weight	34kDa
Shipping	Ice bag