


Anti-Nmnat3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96T66

Item Number: CM0016738



Introduction

Rabbit polyclonal antibody against human NMNAT3, a mitochondrial enzyme involved in NAD⁺ synthesis, axonal protection, and protein quality control. Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	NMNAT3 Rabbit pAb
Remarks/Alias	PNAT3; FKSG76; PNAT-3; NMNAT3
Species	Human
Gene ID	349565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human NMNAT3 (NP_835471.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96T66
Protein Weight	28kDa
Shipping	Ice bag