

Anti-Nmu Polyclonal Antibody For Wb, Elisa - P48645

Item Number: CM0028402



Introduction

Rabbit polyclonal antibody targeting human NMU (Neuromedin-U) validated for Western blot and ELISA. Shows cross-reactivity with mouse and rat. Suitable for studies on feeding, energy balance, inflammation, and neuropeptide signaling.

[Learn More](#)

Parameter	Value
Product Name	NMU Rabbit pAb
Remarks/Alias	NMU
Species	Human
GeneID (Human)	10874
GeneID	10874
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 35-174 of human NMU (NP_006672.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48645
Protein Weight	20kDa
Shipping	Ice bag