

Anti-Nmur1 Polyclonal Antibody For Wb, Elisa - Q9Hb89

Item Number: CM0020283



Introduction

Anti-NMUR1 rabbit polyclonal antibody validated for Western blot (WB) and ELISA in human and mouse samples. Targeting the G protein-coupled receptor NMUR1 involved in feeding, energy balance, and smooth muscle contraction.

[Learn More](#)

Parameter	Value
Product Name	NMUR1 Rabbit pAb
Remarks/Alias	FM3; FM-3; GPC-R; GPR66; NMU1R; (FM-3); NMUR1
Species	Human
Gene ID	10316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-426 of human NMUR1 (NP_006047.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HB89
Protein Weight	47kDa
Shipping	Ice bag