


Anti-Nnmt Rabbit Monoclonal Antibody For Wb, Elisa - P40261

Item Number: CM0007524



Introduction

Recombinant anti-NNMT rabbit monoclonal antibody validated for western blot and ELISA. Targets human nicotinamide N-methyltransferase (29 kDa) for metabolic and epigenetic research. Ideal for studying methylation and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	NNMT Rabbit mAb
Remarks/Alias	NNMT; nicotinamide N-methyltransferase
Species	Human
Gene ID	4837 (human)
Immunogen	Recombinant fusion protein corresponding to amino acids 1-264 of human NNMT (NP_006160.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P40261
Protein Weight	29 kDa
Shipping	Ice bag