


Anti-Nnmt Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P40261

Item Number: CM0004088



Introduction

Anti-NNMT Rabbit Monoclonal Antibody, validated for KD, reacts with human NNMT, suitable for WB, IF/ICC, IP, ELISA. Useful for studying nicotinamide N-methyltransferase function in metabolism and epigenetics.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] NNMT Rabbit mAb
Remarks/Alias	Nicotinamide N-methyltransferase
Species	Human
GeneID (Human)	4837
GeneID	4837
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human NNMT (NP_006160.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P40261
Protein Weight	30kDa
Shipping	Ice bag