


Anti-Nnmt Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P40261

Item Number: CM0016957



Introduction

Rabbit polyclonal antibody against human NNMT, validated for WB, IF-P, ELISA. Reacts with human, mouse, and rat. For studying nicotinamide metabolism, epigenetic regulation, and metabolic diseases.

[Learn More](#)

Parameter	Value
Product Name	NNMT Rabbit pAb
Remarks/Alias	NNMT
Species	Human
GeneID (Human)	4837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human NNMT (NP_006160.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40261
Protein Weight	30kDa
Shipping	Ice bag