


Anti-Nnmt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P40261

Item Number: CM0029593



Introduction

Rabbit polyclonal antibody against human nicotinamide N-methyltransferase (NNMT, 30 kDa). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for epigenetic, metabolic, and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	NNMT Rabbit pAb
Remarks/Alias	NNMT
Species	Human
GeneID (Human)	4837
GeneID	4837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human NNMT (NP_006160.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40261
Protein Weight	30kDa
Shipping	Ice bag