


Anti-Nnmt Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P40261

Item Number: CM0021545



Introduction

High-quality rabbit polyclonal antibody targeting nicotinamide N-methyltransferase (NNMT). Validated for WB, IHC-P, IF-P, and ELISA in human, mouse, and rat samples. Supports research in methylation, epigenetics, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	NNMT Rabbit pAb
Remarks/Alias	NNMT
Species	Human
GeneID (human)	4837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human NNMT (NP_006160.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40261
Protein Weight	30kDa
Shipping	Ice bag