


Anti-Nnmt Rabbit Polymab® For Wb - P40261

Item Number: CM0006999



Introduction

Recombinant rabbit monoclonal antibody against human NNMT, suitable for Western blot, IF/ICC, immunoprecipitation, and ELISA. Cross-reacts with mouse and rat. Useful for research on nicotinamide metabolism and epigenetic regulation.

[Learn More](#)

Parameter	Value
Product Name	NNMT Rabbit PolymAb®
Remarks/Alias	Nicotinamide N-methyltransferase
Species	Human
GeneID (Human)	4837
GeneID	4837
Immunogen	Recombinant fusion protein corresponding to amino acids 1-264 of human NNMT (NP_006160.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40261
Protein Weight	30kDa
Shipping	Ice bag