

Anti-Nnos Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P29475

Item Number: CM0010913



Introduction

Rabbit monoclonal antibody against human nNOS (Nitric oxide synthase 1, P29475). Validated for WB, IF/ICC, IF-P, ELISA with cross-reactivity to mouse and rat. Suitable for neuroscience and NO signaling research.

[Learn More](#)

Parameter	Value
Product Name	nNOS Rabbit mAb
Remarks/Alias	NOS; bNOS; nNOS; IHPS1; N-NOS; NC-NOS
Species	Human
GeneID (Human)	4842
GeneID	4842
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1335-1434 of human nNOS (P29475).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29475
Protein Weight	161kDa
Shipping	Ice bag