


## Anti-Nob1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9ULX3

Item Number: CM0024576



### Introduction

Affinity-purified rabbit polyclonal antibody against human NOB1, validated for Western blot (WB), IF/ICC, and ELISA. Recognizes human, mouse, and rat NOB1 (UniProt Q9ULX3) for ribosome biogenesis research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | NOB1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | ART-4; NOB1P; MST158; MSTP158; PSMD8BP1; NOB1                                                                   |
| Species          | Human                                                                                                           |
| GeneID           | 28987                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human NOB1 (NP_054781.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | Q9ULX3                                                                                                          |
| Protein Weight   | 47kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |