

Anti-Noc4L Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bvi4

Item Number: CM0014601



Introduction

Rabbit polyclonal antibody targeting human NOC4L, validated for Western blot and ELISA. Detects ~58 kDa nucleolar complex protein 4 homolog localized to nucleolus and nucleus membrane. Ideal for nucleolar biology studies.

[Learn More](#)

Parameter	Value
Product Name	NOC4L Rabbit pAb
Remarks/Alias	NOC4; NET49; UTP19; NOC4L
Species	Human
Gene ID (Human)	79050
Gene ID	79050
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human NOC4L (NP_076983.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BVI4
Protein Weight	58kDa
Shipping	Ice bag