

Anti-Nod2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Hc29

Item Number: CM0019551

Introduction

High-quality anti-NOD2 rabbit polyclonal antibody validated for Western blot, immunofluorescence/ICC, and ELISA. Targets human NOD2 (Q9HC29) involved in innate immunity and inflammatory bowel disease.

[Learn More](#)

Parameter	Value
Product Name	NOD2 Rabbit pAb
Remarks/Alias	CD; ACUG; BLAU; IBD1; YAOS; BLAUS; NLRC2; NOD2B; CARD15; CLR16.3; PSORAS1; NOD2
Species	Human
GeneID (Human)	64127
GeneID	64127
Immunogen	Recombinant protein corresponding to a sequence within amino acids 28-270 of human NOD2 (NP_071445.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HC29
Protein Weight	115kDa
Shipping	Ice bag