

Anti-Nod2 Rabbit Polyclonal Antibody For Wb - Q9Hc29

Item Number: CM0024653

Introduction

Rabbit polyclonal antibody against human NOD2 (CARD15). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Detects a band at 115kDa. Suitable for innate immunity and inflammatory bowel disease research.

[Learn More](#)

Parameter	Value
Product Name	NOD2 Rabbit pAb
Remarks/Alias	CD; ACUG; BLAU; IBD1; YAOS; BLAUS; NLRC2; NOD2B; CARD15; CLR16.3; PSORAS1; NOD2
Species	Human
GeneID (human)	64127
GeneID	64127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-216 of human NOD2 (NP_071445.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC29
Protein Weight	115kDa
Shipping	Ice bag