

Anti-Nodal Polyclonal Antibody For Wb, Elisa - Q96S42

Item Number: CM0028535



Introduction

Rabbit polyclonal antibody against NODAL, validated for WB and ELISA. Reacts with human and rat. Essential for embryonic development mesoderm formation. Ideal for developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	NODAL Rabbit pAb
Remarks/Alias	HTX5; NODAL
Species	Human
Gene ID (Human)	4838
Gene ID	4838
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 198-347 of human NODAL (NP_060525.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96S42
Protein Weight	40kDa
Shipping	Ice bag