


Anti-Nogo Receptor 1 (Rtn4R) Monoclonal Antibody For Wb, Elisa - Q9Bzr6

Item Number: CM0010219



Introduction

This Nogo Receptor 1 (RTN4R) rabbit monoclonal antibody is validated for WB and ELISA, recognizing human and mouse targets. Ideal for axon regeneration and neuronal plasticity research.

[Learn More](#)

Parameter	Specification
Product Name	Nogo Receptor 1 (RTN4R) Rabbit mAb
Remarks/Alias	NGR; NOGOR; Nogo Receptor 1 (RTN4R)
Species	Human
Gene ID (Human)	65078
Gene ID	65078
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human Nogo Receptor 1 (RTN4R) (Q9BZR6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZR6
Protein Weight	51kDa
Shipping	Ice bag