


Anti-Nogo Receptor 1 (Rtn4R) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bzr6

Item Number: CM0027598



Introduction

Rabbit polyclonal antibody against human NOGO Receptor 1 (RTN4R), validated for WB, IF/ICC, ELISA. Detects Human, Mouse, Rat. Suitable for axon regeneration and neuronal plasticity studies.

[Learn More](#)

Parameter	Specification
Product Name	NOGO Receptor 1 (RTN4R) Rabbit pAb
Remarks/Alias	NGR; NOGOR; NOGO Receptor 1 (RTN4R)
Species	Human
Gene ID (Human)	65078
Gene ID	65078
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 244-473 of human NOGO Receptor 1 (RTN4R) (NP_075380.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZR6
Protein Weight	51kDa
Shipping	Ice bag