


Anti-Nogo (Rtn4) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Nqc3

Item Number: CM0007264



Introduction

High-quality rabbit monoclonal antibody targeting human Nogo (RTN4, Q9NQC3). Validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for ER and neuroscience research, providing reliable performance in diverse applications.

[Learn More](#)

Parameter	Value
Product Name	Nogo Rabbit mAb
Remarks/Alias	ASY; NSP; NOGO; RTN-X; NSP-CL; RTN4-A; RTN4-C; RTN4-B1; RTN4-B2; NI220/250; Nbla00271; Nbla10545; Nogo
Species	Human
Gene ID (Human)	57142
Gene ID	57142
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 60-160 of human Nogo (NP_065393.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQC3
Protein Weight	130kDa
Shipping	Ice bag