

Anti-Nol12 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ugy1

Item Number: CM0013595

Introduction

NOL12 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Recognizes human NOL12, a nucleolar RNA-binding protein involved in DNA stress resolution and ribosome biogenesis. Suitable for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	NOL12 Rabbit pAb
Remarks/Alias	Nop25; dj37E16.7; NOL12
Species	Human
GeneID (Human)	79159
Gene ID	79159
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-213 of human NOL12 (NP_077289.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UGY1
Protein Weight	25kDa
Shipping	Ice bag