


## Anti-Nol3 Polyclonal Antibody For Wb, If/Icc, Elisa - O60936

Item Number: CM0015305



### Introduction

NOL3 Rabbit pAb detects human and mouse NOL3, a 23 kDa apoptosis repressor. Validated for WB, IF/ICC, ELISA. Suitable for apoptosis and cell death research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | NOL3 Rabbit pAb                                                                                                 |
| Remarks/Alias    | ARC; FCM; MYP; NOP; NOP30; MYOCL1; NOL3                                                                         |
| Species          | Human                                                                                                           |
| Gene ID (Human)  | 8996                                                                                                            |
| Gene ID          | 8996                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human NOL3 (NP_003937.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse                                                                                                    |
| SWISS            | O60936                                                                                                          |
| Protein Weight   | 23kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |