


Anti-Nol3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60936

Item Number: CM0028060



Introduction

Rabbit polyclonal antibody against human Nucleolar protein 3 (NOL3/ARC). Validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Ideal for apoptosis and cell death studies.

[Learn More](#)

Parameter	Value
Product Name	NOL3 Rabbit pAb
Remarks/Alias	ARC; FCM; MYP; NOP; NOP30; MYOCL1; NOL3
Species	Human
GeneID (Human)	8996
GeneID	8996
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human NOL3 (NP_003937.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60936
Protein Weight	23kDa
Shipping	Ice bag