


Anti-Nol6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H6R4

Item Number: CM0018661



Introduction

Rabbit polyclonal NOL6 antibody validated for WB and ELISA. Recognizes human, mouse, and rat NOL6, a 128 kDa nucleolar protein involved in ribosome biogenesis. Ideal for SSU processome studies.

[Learn More](#)

Parameter	Value
Product Name	NOL6 Rabbit pAb
Remarks/Alias	NRAP; UTP22; bA311H10.1; NOL6
Species	Human
GeneID	65083
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1032-1146 of human NOL6 (NP_075068.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6R4
Protein Weight	128kDa
Shipping	Ice bag