

Anti-Nolc1 Rabbit Polyclonal Antibody For Wb - Q14978

Item Number: CM0029212

Introduction

Rabbit polyclonal antibody targeting human NOLC1 (Nopp140). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for nucleolar function and ribosome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	NOLC1 Rabbit pAb
Remarks/Alias	P130; Srp40; NOPP130; NOPP140; NSSATP13; NOLC1
Species	Human
GeneID (Human)	9221
GeneID	9221
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-699 of human NOLC1 (NP_004732.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14978
Protein Weight	74kDa
Shipping	Ice bag