

Anti-Nomo1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15155

Item Number: CM0016266

Introduction

Rabbit polyclonal antibody targeting Nomo1, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for endoplasmic reticulum multi-pass translocon and membrane protein insertion research.

[Learn More](#)

Parameter	Value
Product Name	NOMO1 Rabbit pAb
Remarks/Alias	PM5; Nomo; NOMO1
Species	Human
Gene ID (Human)	23420
Gene ID	23420
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-1000 of human NOMO1 (NP_055102.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15155
Protein Weight	134kDa
Shipping	Ice bag