

Anti-Nomo2 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Jpe7

Item Number: CM0020415



Introduction

Anti-NOMO2 Rabbit Polyclonal Antibody (Nodal modulator 2, pM5 protein 2) for Western Blot (WB) and ELISA. Detects human NOMO2 (139 kDa), a component of the multi-pass translocon complex. Immunogen: recombinant protein (aa 873-1222).

[Learn More](#)

Parameter	Value
Product Name	NOMO2 Rabbit pAb
Remarks/Alias	PM5; Nomo; NOMO2
Species	Human
GeneID (human)	283820
GeneID	283820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 873-1222 of human NOMO2 (NP_001004060.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5JPE7
Protein Weight	139kDa
Shipping	Ice bag